

FILTERING FACEPIECE APPROVAL LABEL



FIDO MASKS CO., LTD. P.O.BOX 101-70, Taipei, Taiwan, R. O. C.

THIS RESPIRATOR IS APPROVED ONLY IN THE FOLLOWING CONFIGURATION

T.C.	Protection	Respirator	Cautions and Limitations
		F720C	
84A-4228	N95	X	ABCJMNOP

1. Protection

N95-Particulate Filter (95% filter efficiency level) effective against particulate aerosols free of oil; time use restrictions may apply

2. Cautions and Limitations

A-Not for use in atmospheres containing less than 19.5% oxygen.

B-Not for use in atmospheres immediately dangerous to life or health.

C-Do not exceed maximum use concentrations established by regulatory standards.

J-Failure to properly use and maintain this product could result in injury or death.

M-All approved respirators shall be selected. Fitted, used, and maintained In accordance with MSHA, OSHA, and other applicable regulations.

N-Never substitute, modify, add, or omit parts. Use only replacement parts In the configuration as specified by the manufacturer.

O-Refer to users instructions, and/or maintenance manuals for information On use and maintenance of these respirators.

P-NIOSH does not evaluate respirators for use as surgical masks.



NIOSH Reference: TN-14223 Mfr. Reference: FISF720C Centers for Disease Control and Prevention (CDC) National Institute for Occupational Safety and Health (NIOSH) National Personal Protective Technology Laboratory (NPPTL)

Pittsburgh, PA 15236-0070 Phone: 412-386-4000 Fax: 412-386-4051 December 21, 2005

P.O. Box 18070

Mr. Arnold Chang Fido Masks Co., Ltd. P.O. Box 101-70 Taipei, TAIWAN, R.O.C.

Dear Mr. Chang:

The National Institute for Occupational Safety and Health (NIOSH) has reviewed your request dated September 24, 2005 and accepted October 6, 2005. This request was for an extension of approval to TC-84A-4228 to add the model F720C N95 filtering facepiece that is similar in design to the model F710C but incorporates plastic hangers for the head strap attachment method.

This request is granted. Approvals are granted only for documentation written in the English language. It is the manufacturer's responsibility to correctly translate materials desired in languages other than English.

The diskette enclosed with this letter contains the final respirator approval label. The abbreviated label has been accepted as submitted. The cautions and limitations which apply to this approval are on the approval label. Only those assemblies affected by this request, or where new approval numbers are assigned, apply to this approval action. Production approval labels cannot include information on unapproved configurations.

The approved assemblies consist of the parts as listed on the approval label and the assembly matrix. Parts are to be marked with the numbers indicated on the approval label in a legible and permanent manner (marking cannot be removed without evidence of its previous presence).

No changes may be made to any respirators and accompanying documentation without prior written approval of NIOSH. Requests for changes must be submitted to NIOSH and a modification of this approval must be granted before changes are made.

Page 2 - Mr. Arnold Chang - TN-14223

Your check, number T0003040477 in the amount of \$100.00 is also enclosed.

Sincerely,

Heinz Ahlers Acting Chief

Technology Evaluation Branch

National Personal Protective Technology Laboratory

Enclosures

National Institute for Occupational Safety and Health Technology Evaluation Branch Certification, Evaluation and Testing Section

Project Background Report

TN-14223

Mfr.: Fido Masks Co., LTD.

Prepared by: Tom Pouchot

Report Date: December 16, 2005

I. Background

In an application dated September 24, 2005 Fido Masks Co. requested an extension of approval to TC-84A-4228 to add the model F720C N95 filtering facepiece that is similar in design to the model F710C but incorporates plastic hangers for the head strap attachment method.

This F710C filtering facepiece model had been previously tested and approved successfully under TN-13807.