



**FILTERING FACEPIECE
APPROVAL LABEL**



FIDO MASKS CO., LTD.
P.O.BOX 101-70, Taipei,
Taiwan, R. O. C.

THIS RESPIRATOR IS APPROVED ONLY IN THE FOLLOWING CONFIGURATION

T.C.	Protection	Respirator	Cautions and Limitations
		REX-350V	
84A-4230	N95	X	ABCJMNOP

1. Protection

N95-Particulate Filter (95% filter efficiency level) effective against particulate aerosols free of oil; time use restrictions may apply

2. Cautions and Limitations

A-Not for use in atmospheres containing less than 19.5% oxygen.

B-Not for use in atmospheres immediately dangerous to life or health.

C-Do not exceed maximum use concentrations established by regulatory standards.

J-Failure to properly use and maintain this product could result in injury or death.

M-All approved respirators shall be selected. Fitted, used, and maintained

In accordance with MSHA, OSHA, and other applicable regulations.

N-Never substitute, modify, add, or omit parts. Use only replacement parts

In the configuration as specified by the manufacturer.

O-Refer to users instructions, and/or maintenance manuals for information

On use and maintenance of these respirators.

P-NIOSH does not evaluate respirators for use as surgical masks.